

Multiplex IHC & Pathology applications

Anastasia Logothetidou
Senior Assay Development Scientist
HistoPathology

25-SEP-23

PATHOLOGY SYMPOSIUM



CellCarta

MAPPING PRECISION MEDICINE

AGENDA

1. CellCarta Histopathology
2. Concept & principle of multiplexing
3. CD8 multiplexing examples
 - CD8/panCK
 - CD8/GZMB/FOXP3

CellCarta: global CRO for end-to-end precision medicine



Clinical Operations














Assay Development



RnD



Glimpse into the multiplex portfolio (IHC, mIHC, mIF, RNA ISH, FISH)




CONCEPT OF MULTIPLEXING

Doing more with less:

Simultaneous detection of multiple markers while preserving precious samples


- Spatial localization (identification & interpretation within the morphology context)
- cellular function/activation and
- cell-cell interactions

PRINCIPLE OF MULTIPLEXING



Use of automated staining platforms and detection kits lead to straightforward & robust assays


- Sequential staining provides flexibility to combine markers
 - Automated denaturing prevents cross reactivity
- Open or ready-to-use detection allows assay customization





CD8/panCK

Numerous studies have described the relationship between the presence and spatial distribution of TILs and treatment benefit

- highlights the role of CTL biomarkers




CD8/panCK



CD8+ lymphocytes stain brown with DAB

Tumor cell cyokeratin stains purple with DISC chromogen


Hematoxylin counterstain



SCORING

Image analysis

Single Section
Topological Analysis
(SSTA)



Pathologist scoring

Density Proportion
Score
(DPS)

IMAGE ANALYSIS WORKFLOW

H&E
EVALUATION
BY
PATHOLOGIST







ANNOTATION
OF TUMOR AND
INVASIVE
MARGIN ON
H&E IMAGE BY
PATHOLOGIST




DETAILED
ANNOTATION
OF ROIS ON
IHC IMAGE
(EXCLUSION
OF ARTIFACTS)




RUNNING TEST
SPECIFIC
ANALYSIS &
QUALITY
CHECK




PanCK - image



PanCK-CD8 – mask



PanCK-CD8 – CD8 detection



Reportables:

- CD8 area relative to
- CD8 density within
- ...

Central tumor
Invasive Margin
Tumor stroma
Carcinoma cell compartment

PATHOLOGIST (DENSITY PROPORTION SCORING) WORKFLOW

SURVEY SLIDE
AT 10X



TUMOR
STROMA (ITS)



TUMOR
NESTS (IE)



PHENOTYPE






- Determine the general distribution of tumor cell nests, cell strands and stroma


- Score the % Immune Cells (IC) in each density bin (absent, low, moderate, high)

- Score the % IC in each density bin (absent, low, moderate, high)

- Combine scores to determine the overall immune phenotype (Desert – Excluded – Inflamed)

DENSITY PROPORTION SCORE







CD8-GZMB-FOXP3

By combining CD8 with the cytolytic GZMB, the density and proportion of the activated CTLs can be assessed



FOXP3+ Tregs are immune suppressive cells, and their abundance/distribution in tumor infiltrates can be associated with unfavorable clinical prognosis



EPITHELIAL NESTS (LOW VS HIGH)



STROMA (LOW VS HIGH)



Phenotype	Image
Desert	
Excluded	
Inflamed	

CD8/GZMB/FoxP3


Excluded

Inflamed

Phenotype = Inflamed (25%) + excluded → inflamed

APPLICATIONS

The IC phenotype is dynamic and will change as a tumor develops and in response to therapy. It is, therefore, a useful parameter to support treatment strategies at different stages of carcinogenesis.



CONCLUDING THOUGHTS

- As the transition from chemotherapy to targeted and immune therapies continues to develop, investigators will utilize more complex clinical trial designs to predict response or resistance to therapy regimens
- Multi-parametric approaches are needed to identify patients most likely to respond to ICI treatment
 - combination of targeted and immune therapies
 - combination of different ICI drugs
- Focus on TME characterization and analysis of the interaction between tumor and effector IC and between effector cells and APCs
 - Spatial relationships and distributions
 - Expression profiles (suggest activation or suppression/exhaustion of effector cells or may hint at mechanisms of resistance)
- **Important to develop robust tissue-based multiplex technologies that can be adopted for routine clinical practice** and – if necessary – successfully shepherded through regulatory approval.